

1. Experiment descriptions and data format

Two-channel PPG signals, three-axis acceleration signals, and one-channel ECG signals were simultaneously recorded from subjects with age from 18 to 35. For each subject, the PPG signals were recorded from wrist by two pulse oximeters with green LEDs (wavelength: 515nm). Their distance (from center to center) was 2 cm. The acceleration signal was also recorded from wrist by a three-axis accelerometer. Both the pulse oximeter and the accelerometer were embedded in a wristband, which was comfortably worn. The ECG signal was recorded simultaneously from the chest using wet ECG sensors. All signals were sampled at 125 Hz and sent to a nearby computer via Bluetooth.

Each dataset with the similar name 'DATA_01_TYPE01' contains a variable 'sig'. It has 6 rows. The first row is a simultaneous recording of ECG, which is recorded from the chest of each subject. The second row and the third row are two channels of PPG, which are recorded from the wrist of each subject. The last three rows are simultaneous recordings of acceleration data (in x-, y-, and z-axis).

During data recording, each subject ran on a treadmill with changing speeds. For datasets with names containing '**TYPE01**', the running speeds changed as follows:

rest(30s) -> 8km/h(1min) -> 15km/h(1min) -> 8km/h(1min) -> 15km/h(1min)
-> rest(30s)

For datasets with names containing '**TYPE02**', the running speeds changed as follows:

rest(30s) -> 6km/h(1min) -> 12km/h(1min) -> 6km/h(1min) -> 12km/h(1min)
-> rest(30s)

For each dataset with the similar name 'DATA_01_TYPE01', the ground-truth of heart rate can be calculated from the simultaneously recorded ECG signal (i.e. the first row of the variable 'sig'). For convenience, we also provide the calculated ground-truth heart rate, stored in the datasets with the corresponding name, say 'DATA_01_TYPE01_BPMtrace'. In each of this kind of datasets, there is a variable 'BPM0', which gives the BPM value in every 8-second time window. Note that two successive time windows overlap by 6 seconds. Thus the first value in 'BPM0' gives the calculated heart rate ground-truth in the first 8 seconds, while the second value in 'BPM0' gives the calculated heart rate ground-truth from the 3rd second to the 10th second.

(NOTE: Data may have noise because of excessive motion during exercise and loose contact between skin and wearable device.)

2. Testing Data Input Format

Each dataset contains a variable named 'sig' of 5 rows. The first and the second rows are the two-channel PPG signals, while the other rows are the tri-acceleration signals.

(NOTE: In the testing datasets, ECG signals are not present)

Signal Sampling Rate: 125 Hz.

3. Estimation Requirement (Output Format)

An algorithm should give the estimated heart rate in each time window of 8 seconds. Two successive time windows overlap by 6 seconds. Below are few examples.

- The first estimated heart rate value should be obtained from the data with the sampling indexes from the 1st to the 1000th.
- The second estimated heart rate value should be obtained from the data with the sampling indexes from the 251th to the 1250th. However, you can also utilize information from the previous time window (such as the previously estimated heart rate). **But note that you cannot use any data with the sampling indexes > 1250th (or associated information).**